

MORADABAD INSTITUTE OF TECHNOLOGY
MORADABAD

Student Name :- Aryabh Prayabati

Father's Name :- Mr. Sanjay Baby

Subject :- Chemistry Practical (KAS-102)

Branch :- C.S.





Course :- B.Tech. First year

Session :- 2019-20

Subject Teacher :- Dr. Harendra Kumar Sharma



I n d e x

S. No.	Name of the Experiment	Page No.	Date of Experiment	Date of Submission	Remarks
1-	To determine the alkalinity in water sample	19	10/8/19	26/8/19	(9)  26/8/19
2-	To determine the temporary and permanent hardness of water by using complexometric titration.	26	10/8/19	02/09/19	(9)  02/09/19
3-	To determine iron(ferrous) content in the supplied Sample of iron ore by titrimetric analysis against Standard $K_2Cr_2O_7$ solution $K_3[Fe(CN)_6]$ as an external indicator	21	09/19	9/9/19	(8)  9/9/19
4-	To determine the viscosity of a liquid with the help of Ostwald's Viscometer	9	10/9/19	16/9/19	(B)  16/9/19
5	To study variation of surface tension of liquid on addition of soap or surfactant of known conc.	16	10/9/19	23/9/19	(B)  23/9/19
6-	To determine chloride content in given water sample by Mohr's method	23	10/9/19	30/9/19	(9)  30/9/19

I n d e x

S. No.	Name of the Experiment	Page No.	Date of Experiment	Date of Submission	Remarks
7	To determine the percentage of available Chlorine in a given sample of bleaching powder by iodometric titration method		30/09/19	(9)	(P) 14/10/19
8	To prepare phenol formaldehyde resin Bakelite. To prepare urea formaldehyde (UF) Resin.		14/10/19	(9)	(P) 4/11/19
9-	To detect the elements and func. Groups in the given organic compound.		4/11/19	(9)	(P) 21/10/19
10-	To determine the Strength of given organic comp' ^{acid} hydrochloric soln on by titrating it against a pH metrically against a Sodium hydroxide		21/10/19	(8)	(P) 21/10/19

General instructions

- ① Laboratory coat must be worn at all times during the laboratory period.
- ② Do not perform unauthorized experiment.
- ③ Never work alone in the laboratory.
- ④ At once report all incidents to your teacher.
- ⑤ Never taste a chemical or soln unless specifically directed to do so.
- ⑥ Be alert to what your laboratory neighbour is doing. Your neighbour accident may be your injury.
- ⑦ Read the label carefully before removing a chemical from its shelf and keep the bottle at its appropriate place.
- ⑧ If any of your glass equipment is cracked or broken, do not use it.
- ⑨ If you need to observe the odour of a substance, gently fan a little of the vapour toward with your hand, but keep your face at a safe distance.
- ⑩ When you are heating a test tube, never point its mouth towards yourself or anyone else.
- ⑪ When diluting concentrated acid pour the acid slowly into water with stirring. Never add water to the acid.
- ⑫ Do not casually dispose of chemicals down the drain.
- ⑬ Protect your hands if you have to pick up hot
- ⑭ Do not force rubber stoppers on the glass objects.
- ⑮ Use no open flames near inflammable solvents.

Teacher's Signature _____

- (16) Use fire extinguisher at the time of fire. Before you leave the laboratory make certain that any gas line you used is shut tightly.

Laboratory accidents and first aid :-


- ⇒ Fire - ① For large fires get yourself and others away from burning area.
 ② Notify your instructor immediately.
 ③ Remove container or flammable solvent from the immediate area and turn off burners.
 ④ If your clothing on fire, do not run, roll on the floor to smother the fire and to help keep flames away from your hands.

⇒ Small Burns : ~~Immediately plunged the burned area into cold water preferably with ice in it.~~

⇒ Severe Burns : Do not apply any ointment but seek professional medical treatment at once.

- ⇒ Chemical Burns :-
- ① General Burns :- Immediately wash off it with plenty of water.
 - ② Acids = Wash with plenty of water then with very dilute sodium bisulphite solution.
 - ③ Bases = Wash with plenty of water then with very dilute acetic acid soln.

Teacher's Signature _____



Experiment - I

- Object = To determine the amount of alkalinity of supplied water sample.
- Apparatus used = ~~Burette~~, pipette, conical flask, funnel, Beaker (250 ml).
- Chemicals used = N/20 HCl, phenolphthalein & methyl orange.
- Theory - "The alkalinity of water may be defined as its capacity to neutralize acids". The alkalinity of water is mainly due to presence of hydroxide ions (OH^-), carbonate ions (CO_3^{2-}) & bicarbonate ions (HCO_3^-) in water.

These can be obtained separately by the titration against standard acid using phenolphthalein & methyl orange as indicators. The various rxn involved are:

The titration of water samples against standard acid up-to phenolphthalein end point marks the completion of rxn (i) and (ii) only. The volume of acids (A ml) used up-to this point corresponds to the complete neutralization of OH^- ions and one half of the normal CO_3^{2-} ions present. The alkalinity up-to phenolphthalein end point is known as phenolphthalein alkalinity (P).

Teacher's Signature _____

The titration of water sample against the standard acid up-to methyl orange end point shows the completion of rxns (i) (ii) and (iii). ∵ the amount of acid (B ml) carbonate and all the bicarbonates present. The total amount of acid $[A+B]$ ml in two titrations represents the total alkalinity (M)

→ The possible combinations of ions causing alkalinity in water are :-

- 1- OH^- only
- 2- CO_3^{2-} only
- 3- HCO_3^- only
- 4- OH^- and CO_3^{2-} only
- 5- CO_3^{2-} and HCO_3^- only

The possibility of OH^- and HCO_3^- together is not possible since they combine together to form CO_3^{2-} and H_2O



After knowing the values of P and M the nature of alkalinity can be determined by using the table

Results of titration to phenolphthalein
end point [P] & Methyl orange end point [M]

$$[P] = 0$$

$$[P] = [M]$$

	OH^-	CO_3^{2-}	HCO_3^-
	Nil	Nil	[M]
	[P] or [M]	Nil	Nil.

Teacher's Signature _____

$[P] = 1/2 [M]$	Nil	$2[P] \text{ or } [M]$	Nil
$[P] > 1/2 [M]$	$2[P] - [M]$	$2[M] - [P]$	Nil
$[P] < 1/2 [M]$	Nil	$2[P]$	$[M] - 2[P]$

→ Alkalinity is generally expressed as parts per million (ppm) in terms of CaCO_3 .

- Observations :-

S.no	Volume of water Sample (ml)	Vol. of acid used upto phenolphthalein End point [P]	Vol. of acid used upto methyl orange End point [M] [ml]
1 -	25 ml	3.9 ml	9.5 ml
2 -	25 ml	3.6 ml	8.0 ml
3 -	25 ml	3.5 ml	8 ml

Calculations :-

(i) Phenolphthalein alkalinity $[P]$ in terms of CaCO_3 equivalent

$$\frac{N_1 V_1}{(Acid)} = \frac{N_2 V_2}{(Water Sample)}$$

$$\frac{1 \times 3.6}{20} = N_2 \times 25$$

$$N = \frac{1 \times 3.6}{500}$$

Strength in terms of CaCO_3 equivalent
 $= N_2 \times \text{Equivalent weight of } \text{CaCO}_3$

Teacher's Signature _____

$$= \frac{3.6}{500} \text{ gm/l}$$

$$= .36 \text{ gm/l}$$

$$= .36 \times 1000 \text{ mg/l or ppm}$$

$$= \underline{360 \text{ mg/l or ppm}}$$

(ii) Methyl orange alkalinity [M] in terms of CaCO_3 equivalent.

$$\text{N}_3 \text{V}_3 = \text{N}_4 \text{V}_4 \\ (\text{Acid}) \quad (\text{Water Sample})$$

$$\frac{1}{20} \times 8 = \text{N}_4 \times \text{V}_4$$

$$\cancel{\frac{1}{20} \times 8} = \text{N}_4 \times 25$$

$$\frac{2}{125} = \text{N}_4$$

Strength in terms of CaCO_3 equivalent -

$$= \text{N}_4 \times \text{Eq. weight of } \text{CaCO}_3$$

$$= \frac{2}{125} \times 50,2$$

$$= \frac{.8}{5} \text{ gm/l}$$

$$= 800 \text{ mg/l or ppm.}$$

- Important Note = To calculate the alkalinity in terms of individual ions, find out the to which bases

Teacher's Signature _____

the values of P and M falls from table A.

- (i) If $P = 0 \rightarrow$ Alkalinity is only due to bicarbonate ions
 \rightarrow Bi-carbonate alkalinity in terms of CaCO_3
 $\text{eq.} = M$
- (ii) If $P = 1/2 M \rightarrow$ Alkalinity is only due to carbonate ions
 \rightarrow Carbonate alkalinity in term of CaCO_3 eq.
 $= 2P$ or M
- (iii) If $P = M \rightarrow$ Alkalinity is only due to hydroxide ions
 \rightarrow Hydroxide alkalinity in terms of CaCO_3
 $\text{eq.} = P$ or M .
- (iv) If $P > 1/2 M \rightarrow$ Alkalinity is due to hydroxide and carbonate ions
a) Hydroxide alkalinity in terms of CaCO_3 eq. = $(2P - M)$
b) Carbonate alkalinity in terms of CaCO_3 eq. = $2(M - P)$
- (v) If $P < 1/2 M \rightarrow$ Alkalinity is due to ~~hydroxide &~~ ^{carbonate} bicarbonate ions
a) Carbonate alkalinity in terms of CaCO_3 eq. = $2P$
b) Bi-carbonate alkalinity in terms of CaCO_3 eq. = $M - 2P$

- Result = The given water sample has :

* alkalinity due to OH^- = Nil

* alkalinity due to CO_3^{2-} = $2 [P] = 2 \times 360 = 720$

$$\text{As } P < \frac{1}{2} M; * \text{ Alkalinity due to } \text{HCO}_3^- = 800 - 2 \times 360 \\ = 80$$

$$\rightarrow \text{Total alkalinity} = \text{Alkalinity due to } (\text{CO}_3^{2-} + \text{HCO}_3^- + \text{OH}^-) \\ = 720 + 80 \\ = 800 \text{ ppm}$$

- Precautions =

- 1) The glass apparatus should be cleaned and properly rinsed before the start of experiment.
- 2) HCl level in burette must be read carefully.
- 3) Volume of HCl should be same in all titrations.

- Viva Voice = Q1- What is alkalinity?

Ans = Alkalinity is the measure of the ability of a solution to neutralize acids.

Q2- Why alkalinity is not due to the simultaneously presence of OH⁻ and CO₃²⁻ ?

Ans - The alkalinity is not due to the OH⁻ and CO₃²⁻ because they combine together to form CO₃²⁻ and H₂O.



Q3- Write the pH range of Phenop. & methyl orange?

Ans - Phenolphthalein = 8.3 - 10.0

Methyl orange = 3.1 - 4.4

9

9/1


Teacher's Signature _____

Experiment - 2

- Object = To determine the temporary or permanent hardness of water by using complexometric titration.
- Apparatus used = Burette, conical flask, funnel, pipette, Beaker, etc.
- Chemicals used = EDTA solⁿ (0.01)M, Buffer Solution ($\text{NH}_4\text{OH} + \text{NH}_4\text{Cl}$), Eriochrome Black-T (EBT), water sample.
- Theory = The hardness of water can be determined by complexometric titration. EDTA is used as complexing agent. The Ca^{2+} and Mg^{2+} present in water are titrated with EDTA using erichrome Black-T as indicator.

Estimation of hardness by EDTA method is based on the following principle.


First, the indicator erichrome Black-T which is blue coloured dye, form an unstable complex with calcium or magnesium ions in hard water at pH of 9 to 10. The complex is wine red in colour.




As this solⁿ is titrated against EDTA, the free Ca^{2+} and Mg^{2+} ions in water form stable metal.

Teacher's Signature _____

ion EDTA complex, which is colourless.



Once the free metal ions are complexed, the EDTA replaced Ca^{2+} and Mg^{2+} ions from the unstable indicator complex also, to form a stable complex with the result, the indicator is set free. Since the free indicator is blue in colour at the above at the above mentioned pH, the end point is the appearance of blue colour.



Thus the amount of EDTA used corresponds to the hardness of water. The temporary hardness of water is removed by boiling and after the removal of ppt by filtration, the permanent hardness in the filtrate is determined by titration with EDTA as above.

Temporary hardness will be given by the difference of total hardness and permanent hardness.

Observations :-

- Total hardness (before boiling)

S.no	Volume of water Sample taken (ml)	Burette Reading			Actual volume of EDTA used
		Initial	Final	Difference	
1	25 ml	0	16	16	
2	25 ml	16	32	16	$a = 16 \text{ ml}$
3	25 ml	32	48.2	16.2	

- Permanent hardness (after boiling)

S.no	Volume of water Sample taken (ml)	Burette Reading			Actual volume of EDTA used.
		Initial	Final	Difference	
1	25 ml	0	14	14	
2	25 ml	14	28	14	$b = 14 \text{ ml}$
3	25 ml	28	40.5	12.5	

- Calculations = EDTA forms 1:1 complex with Ca^{2+} and Mg^{2+}
 $\therefore 1000 \text{ ml of } 1\text{M EDTA} = 1 \text{ mole of } \text{CaCO}_3$
 $= 100 \text{ gm of } \text{CaCO}_3$
 $\therefore 1 \text{ ml of } 0.01 \text{ EDTA} = \frac{100}{1000} \times \frac{1}{100} \text{ gm of } \text{CaCO}_3$

Teacher's Signature _____

$$= 1 \text{ mg of } \text{CaCO}_3$$

$$\text{Total hardness} = 16 \text{ ml of } 0.01 \text{ EDTA} = 1 \times 16$$

$$= 16 \text{ mg of } \text{CaCO}_3$$

This hardness is present in 25 ml of sample water.

So, hardness present in 1000 ml of sample water

$$= 1000 \times \frac{16}{25}$$

$$= 40 \times 16 = 640 \text{ gm of } \text{CaCO}_3$$

$$\therefore \text{Total hardness} = 640 \text{ mg/l lit or ppm}$$

$$\text{Permanent hardness} = 14 \text{ ml of } 0.01 \text{ EDTA} = 1 \times 14 = 14 \text{ mg of } \text{CaCO}_3$$

This hardness is present in 25 ml of sample water

So, hardness present in ~~25 ml~~ 1000 ml of sample water

$$= 1000 \times \frac{14}{25}$$

$$= 40 \times 14 = 560 \text{ gm of } \text{CaCO}_3$$

$$\text{Temporary hardness} = (\text{Total} - \text{Permanent}) \text{ hardness}$$

$$= 640 - 560$$

$$= 80 \text{ mg/l lit or ppm.}$$

• Result Total hardness = 640 ppm

Permanent hardness = 560 ppm

Temporary hardness = 80 ppm

• Precautions

- 1) The glass apparatus should be cleaned before starting the experiment.
- 2) The pH should be maintained during titration.
- 3) Hard water sample should be boiled properly.

• Viva Voice :-

Q1- How will you define hardness?

Ans- The hardness of water is defined as the presence of (Ca^{2+}) and (Mg^{2+}) in water.

Q2- What is the cause of temporary hardness in water?

Ans- Temporary hardness is due to presence of $Ca(HCO_3)_2$ and $Mg(HCO_3)_2$ in water.

Q3- How will you remove the temporary hardness?

Ans- It is removed by heating the water as calcium & magnesium carbonates are precipitated out of water and are collected as residue.

Q4- What do you understand by permanent hardness?

Ans- It is the presence of chlorides, nitrates and sulphates of calcium and magnesium in water.

Q5- What is the effect on the end point with erichrome black-T indicated, if the hard water sample does

Teacher's Signature _____



not contain Mg^{2+} / Ca^{2+} ?

Ans- EBT will not change its colour to wine red.
The soln will be colour blue. Titration will not occur.

(9)

(5) 

Teacher's Signature _____



Experiment - 3

Object = To determine iron (ferrous) content in the supplied sample of iron ore by titrimetric analysis against Standard $K_2Cr_2O_7$ solution using $K_3[Fe(CN)_6]$ as an external indicator.

Apparatus used :- Burette, pipette, conical flask, beaker, glass rod, watch glass, funnel.

Chemicals used = Iron ore sample, standard potassium dichromate solution ($N/10$), potass potassium ferricyanide

Theory = Potassium dichromate solution in presence of H_2SO_4 is a strong oxidizing agent which converts ferrous ions (Fe^{2+}) present in the iron ore into ferric ions (Fe^{3+}). The end point can be determined by using potassium ferricyanide as an external indicator as after the completion of rxn, no Fe^{2+} will be present in the rxn thereby no greenish blue coloration is produced.

Observations :-

Sno	Vol. of Sample Soln	Reading of Burette (ml)		Actual volume of $K_2Cr_2O_7$ used (ml)
		Initial	Final	Difference
1	25 ml	0	12	12
2	25 ml	12	24.5	12.5
3	25 ml	24.5	37	12.5

Teacher's Signature _____

Calculations :-

$$N_1 V_1 = N_2 V_2$$



$$N_1 \times 25 = \frac{N}{10} \times 12.5$$

$$N_1 = \frac{12.5}{25 \times 10}$$

Strength of given iron or soln = Normality \times
Eq. weight of molar salt.

$$= \frac{12.5 \times 392.12}{25 \times 10} \text{ gm/lit.}$$

$$= 19.606 \text{ gm/lit.}$$

As 392.12 gm of Molar salt contain = 56 gm Fe
 So, 19.606 gm of Molar salt contain = $\frac{56 \times 392.12}{19.606}$ gm
 $= 2.8 \text{ gm}$

Result = The ferrous content in the supplied sample of iron ore is 2.8 gm.

- Precautions
- ① Potassium dichromate acts as an oxidising agent in acidic medium, ∴ add dilute H₂SO₄ in molar salt soln.
 - ② The vol. of soln taken out for checking end point should not be more otherwise the result.

Teacher's Signature _____

will not give accurate value.

- (3) The glass rod should be washed each time before withdrawing a drop from the conical flask.
- (4) Rinse the burette and pipette properly before starting the titration.

Viva Voice

Q1- What is the formula of ferrous ammonium sulphate?

Ans- $(\text{NH}_4)_2 \text{Fe}(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$

Q2- What are the main ores of iron?

Ans- Haematite, Magnetite, Limonite, Siderite.

Q3- What types of indicator are using in this titration?

Ans- $\text{K}_3[\text{Fe}(\text{CN})_6]$ as an external indicator.

8

P

Teacher's Signature _____

Experiment - 4

Object = To determine the viscosity of a liquid with the help of a Ostwald's viscometer.

Theory Apparatus Used = Ostwald's viscometer, Specific gravity bottle, stopwatch, beaker, etc.

Theory = Viscosity of liquid is a measure of its frictional resistance. Let 'l' be length of the capillary and 'r' its radius. Then the volume of liquid flowing through the capillary tube for difference of pressure P is given by the formula.

$$V = \frac{\pi P r^4}{8 n l}$$

Where n is the viscosity of a liquid taken in viscometer.

For first liquid, $V_1 = \frac{\pi P_1 r_1^4}{8 n_1 l} \quad \text{--- (1)}$

For second liquid $V_2 = \frac{\pi P_2 r_2^4}{8 n_2 l} \quad \text{--- (2)}$

i.e. $\frac{V_1}{V_2} = \frac{P_1}{P_2} \times \frac{n_2}{n_1} \quad \text{--- (3)}$

Now the pressure of the liquids $P = h \rho g$ or $h \rho g$
Since h and g are same for two liquids.

Teacher's Signature _____

Hence $P \propto \rho$

Pressure of the liquid in each case is proportional to the density of liquid. $P_1 \propto \rho_1$, $P_2 \propto \rho_2$

$$\frac{P_1}{P_2} = \frac{\rho_1}{\rho_2} \quad \textcircled{4}$$

Let Q be the volume of each liquid flowing through the capillary tube in time t_1 and t_2 . Then -

$$V_1 = \frac{Q}{t_1} \quad \& \quad V_2 = \frac{Q}{t_2}$$

$$\frac{V_1}{V_2} = \frac{t_2}{t_1} \quad \textcircled{5}$$

Substituting the value of $\frac{P_1}{P_2}$ and $\frac{V_1}{V_2}$ from $\textcircled{4}$ &

$\textcircled{5}$ respectively in the equation $\textcircled{2}$ we get -

$$\frac{t_2}{t_1} = \frac{\rho_1}{\rho_2} \times \frac{n_2}{n_1}$$

$$\boxed{\frac{n_1}{n_2} = \frac{\rho_1 t_1}{\rho_2 t_2}}$$

Thus the comparison of viscosities is made by knowing the times for corresponding flow and the densities of the liquids.

Teacher's Signature _____

Observations :- @ Table of measurement of flow time
[t_1 and t_2]

S.no	Water		Liquid	
	Flow time b/w A and B	Mean	Flow time b/w A and B	Mean
1	53 sec		65 sec	
2	50 sec	50.2	65 sec	63.8
3	50.8 sec		63 sec	
4	49 sec		64 sec	
5	49 sec		62 sec	

⑥ Table for measurement of density (ρ_2)

Weight of empty RD bottle \rightarrow 21 gm

Weight of RD bottle + water \rightarrow 50.14 gm

Weight of RD bottle + liquid \rightarrow 46.30 gm

Weight of Water filled \rightarrow 29.14

Weight of Liquid filled \rightarrow 25.3


$$\text{Density of liquid } \rho_2 = \frac{25.3}{29.14} = .86$$

⑦ Room temperature $= 20^\circ\text{C}$

⑧ Density of water $= \rho_1 = 1.00$

⑨ Viscosity of water at $20^\circ\text{C} = 0.0101$ poise

Teacher's Signature _____



Ostwald's Viscometer

Calculations :- The viscosity of liquid is calculated from the following formula assuming η_1 and P_1 of water 0.0101 and 1.00 respectively, at $20^\circ C$.

$$\frac{\eta_1}{\eta_2} = \frac{f_1}{f_2} \times \frac{P_1}{P_2}$$

$$\text{viscosity of liquid } \eta_2 = \frac{\eta_1 \times f_2 \times P_2}{t_1 \times P_1}$$

$$= \frac{0.0101 \times 63.8 \times 0.86}{50.2 \times 1.00}$$

$$= \frac{.55}{50.2}$$

$$\eta_2 = 0.0109 \text{ poise.}$$

Result → Viscosity of liquid at $20^\circ C$ is 0.0109 poise.

Precautions →

- ① The flow time of a liquid filled in viscometer should be between 1 to 10 minutes.
- ② The viscometer and RD bottle should be thoroughly cleaned
- ③ The viscometer should be adjusted in any accurately vertical position.

(7) Five or six observations should be recorded in each case.

Viva Voice :-

Q1- What is viscosity ?

Ans- The viscosity of a fluid is a measure of its resistance to flow at a given rate.

Q2- What is unit of viscosity ?

Ans- Poise (dyne s/cm²)

Q3- What is the effect of temperature on viscosity ?

Ans- On high temperature viscosity increases in gases but decreases in liquids.

Q4- Which is more viscous honey or water ?

Ans- Honey because it has more strong intermolecular attractive force than water.

(8)

Teacher's Signature _____

Experiment - 5

Object = To Study Variation of surface tension of liquid on addition of soap or surfactant of known concentration

Apparatus Required = Stalagmometer.

Chemicals Used = (i) Reference liquid (H_2O)
(ii) Experimental liquid (Ethyl alcohol)

Theory = Surface tension = The force in newton acting at right angles along the surface of a liquid one metre in length!

Units = N/m SI system or dynes/cm in CGS system.

Procedure = Observations :-

Density of water ($25^\circ C$) =

Density of liquid ($25^\circ C$) =

Surface tension of water ($25^\circ C$) =

Observation table :-

S.no	Water		Given liquid	
	No. of drops	Mean	No. of drops	Mean
1)	40		72	
2)	36	38	71	71
3)	38		67	
	38		71	

Teacher's Signature _____

Calculations :-

Surface tension of experimental liquid

$$= \frac{\text{No. of drops of H}_2\text{O} \times \frac{4}{7} \times \text{Surface tension}}{\text{No. of drops of exp. liq.}} \text{ of water (dynes/cm)}$$

$$= \frac{38}{71} \times \frac{1.0071089}{1} \times 72.14$$

$$= 38.6915 \text{ dynes/cm}$$

Result :-

The surface tension of given liquid are found to be 38.6915 dyne/cm

Precautions :-

- ① Before use, the stalagmometer should be cleaned and dry.
- ② No air bubble should be formed while sucking the liquid into the stalagmometer.
- ③ Throughout the drop counting the process, the stalagmometer should be held in a vertical position.
- ④ Drop formation should be adjusted at a slower rate and should not exceed fifteen drops per min.

Teacher's Signature _____

Viva-Voce :-

Q2- What is surface energy?

Ans- Potential energy per unit area of the liquid surface under isothermal conditions is known as surface energy.

Q3- What is the effect of temperature on surface tension?

Ans- Surface tension decreases when temp^o increases and vice-versa.

Q4- Give some applications of surface tension?

Ans- A needle placed on water can be made to float due to the surface tension of water.

② Surface tension prevents water from passing through the pores of an umbrella.

③ Antiseptics like dettol have low surface tension so that they spread faster.


Q1- What are the effects of surface tension?

Ans- The effect is an inward force at its surface that causes the liquid to behave as if its surface were covered with a stretched elastic membrane.

(8)

(F)

Teacher's Signature _____



Experiment - 6

Object = To determine the chloride content in given water sample by Mohr's method.

Apparatus Used = Burette, pipette, conical flask, Beaker, volumetric flask.

Chemicals Used = Standard ($N/20$) AgNO_3 soln, Water sample, potassium chromate.

Theory - Chlorides usually occur in water as NaCl , CaCl_2 , MgCl_2 . Chloride content in water sample is determined by Mohr's method. A weakly alkaline method soln (7-8) pH of the water sample is titrated against standard silver nitrate solution using K_2CrO_4 as internal indicator. Cl^- ions present in water sample react with Ag^+ ions (coming from AgNO_3 soln) forming insoluble white PPT of AgCl . CrO_4^{2-} ions from indicator also react with Ag^+ ions giving transient red colour, which immediately disappears due to presence of greater concentration of Cl^- ions.

As the titration proceeds all the chloride ions are removed as AgCl , any extra drop of AgNO_3 reacts with excess potassium chromate giving reddish brown colour of Ag_2CrO_4 . The formation of reddish brown colour indicates the end point.

Teacher's Signature _____

Observations - I Titration with water sample =

S.no	Vol. of water Sample(ml)	Burette Reading (ml)			Vol. of AgNO_3 used
		Initial	Final	Difference	
	25	0	6.5	6.5	
1	25	6.5	13	6.5	
2	25	13	19.5	6.5	6.5 ml
3	25	19.5	26	6.5	
4	25	26			

• Titration with Distilled water (Blank titration)

S.no	Volume of water Sample ml	Burette Reading (ml)			Volume of AgNO_3 used
		Initial	Final	Difference	
1	25	26	27	1.0	
2	25	27	27.5	.5	.5 ml
3	25	27.5	28	.5	
4	25	28	28.8	.8	

Calculations = $N_1 V_1 = N_2 V_2$
 (Water sample) (AgNO_3 soln)

$$\Rightarrow \text{where } V_2 = 6.5 - .5 = 6 \text{ ml}$$

$$\Rightarrow N_1 \times 25 = \frac{1}{20} \times 6 = .012$$

$$\Rightarrow \text{Strength of } \text{Cl}^- \text{ ions} = .012 \text{ } N_1 \times \text{Eq. weight}$$

$$= .012 \times 35.5 \text{ gm/l}$$

$$= .426 \text{ gm/l}$$

$$= 426 \text{ mg/l or ppm}$$

Teacher's Signature _____

Result = The chloride ion conc. in the given water sample is 426 mg/l or ppm

Precautions = ① The glass apparatus should be cleaned with distilled water.

- ② The m/s^n content should be continuously shaken during titration.
- ③ Equal volume of water sample and distilled water should be taken.
- ④ Equal drops of indicator should be used in all titration.
- ⑤ The pH of soln must be maintained b/w 7-8.

Viva Voce

Q1- Name the forms in which Cl^- ions are present?

Ans- NaCl , CaCl_2 , MgCl_2

Q2- Name the indicator used in this titration?

Ans- K_2CrO_4

Q3- What is conc. of Cl^- , which gives unacceptable taste of water?

Ans- Above 250 mg/l conc. of Cl^- it gives salty taste to water.

Q4- What is equivalent weight of Cl^- ion?

Ans- 35.5 gm/l.

Teacher's Signature






Experiment - 7

Object = To determine the percentage of available chlorine in a given sample of bleaching powder.

Chemicals used = Bleaching powder potassium iodide, dilute acid acetic acid, $\text{Na}_2\text{S}_2\text{O}_3$ sodium thiosulphate soln, distilled water and freshly prepared starch soln.

Apparatus = Burette, pipette, measuring flask, conical flask, funnel, pastel and mortar.

Theory = Bleaching powder is famous for its bleaching action. It is also used as disinfectant. The main constituents of bleaching powder are calcium hypochlorite $[\text{Ca}(\text{OCl})_2 \cdot 4\text{H}_2\text{O}]$, the basic calcium chloride $[\text{CaCl}_2 \cdot \text{Ca}(\text{OH})_2 \cdot \text{H}_2\text{O}]$ and some free calcium hydroxide $[\text{Ca}(\text{OH})_2]$. Out of these constituents, the active constituents of bleaching powder is hypochlorite. Hypochlorite may be represented as:



When dilute acid reacts with bleaching powder then free chlorine is liberated acc. to the following rxn-

The amount of chlorine that can be obtained by action of dilute acids on bleaching powder is known as available chlorine. It is expressed as %.

Teacher's Signature _____

weight of bleaching powder.

Observations :

- (a) Weight of empty weight tube = 3.190 gm
- (b) Weight of weighing tube + weight of bleaching powder = 6.75 gm
- (c) Weight of bleaching powder = 3.56 gm

Sno	Vol. of Bleaching Powder.	Burette Initial	Reading (ml)	Difference	Vol. of Hypo used (ml)
1	25 ml	0	38.5	38.5	
2	25 ml	0	38	38	38 ml
3	25 ml	0	37.5	37.5	
4	25 ml	0	38	38	

Calculations

$$N_1 V_1 = N_2 V_2$$

(Bleaching powder) (Sodium thiosulphate)

$$\text{Normality of available chlorine} = N_1 = (N_2 \times V_2) / V_1$$

$$N_1 = \frac{N \times 38}{25} = 0.152$$

$$z = \text{Amount of available chlorine} = N_1 \times 38.5 = 3.55 \times 0.152 \\ \text{per } 250 \text{ ml of soln} = 5.396 \text{ gm/250 ml}$$

$$\text{Amount of avail. chlorine/litre of soln} = \frac{5.396 \times 1000}{250} = 21.584 \text{ gm/lit} \\ = 1.349 \text{ gm/lit}$$

$$\text{Percentage of available chlorine} = \frac{z}{w} \times 100 = \frac{5.396}{3.56} \times 100 = 151.47\% \\ [w = (\text{weight of bleaching powder})] \text{ Teacher's Signature } \underline{\hspace{10cm}}$$

Result = The percentage of available chlorine present in given sample of bleaching powder sample is 37.47%

Precautions - (1) Sample of bleaching powder must not be kept as open as it absorb moisture from the atmosphere.

(2) No lumps of bleaching powder should remain while making the paste.

(3) Rinse the burette with standard $\text{N}/10$ sodium thiosulphate soln.

(4) The starch soln and acetic acid must be added in equal amount in each experiment.

(5) Titration must be started immediately after the addition of acid because chlorine is liberated instantly.

Viva Voce

Q1- What do you mean by available chlorine?

Ans- The amount of chlorine that can be obtained by the action of dilute acids on bleaching powder is known as available chlorine.

Q2- What do you understand by iodometry?

Ans- It is a iodometric titration, a redox titration where the appearance or disappearance of elementary.

Teacher's Signature _____

iodine indicates the end point.

Q3- What is the other name of bleaching powder?

Ans- Calcium Hypochloride

Q4- What is the app. percentage of available chlorine in commercial bleaching powder?

Ans- ~~35%~~ 10 - 40%

Q5- Give the names of main constituents of bleaching powder?

Ans- Calcium hypochloride chloride [Ca(OCl) Cl. 4H₂O]
 Basic calcium chloride [CaO₂. Ca(OH)₂. H₂O]
 free calcium hydioxide [Ca(OH)₂]




⑨ ~~Par~~

Teacher's Signature _____

Experiment - 8

Object - To prepare phenol formaldehyde resin bakelite.

Requirements - Beaker, measuring cylinder, weight box, chemical balance, Filter paper, Formaldehyde, Glacial acetic acid, Phenol, conc. HCl

Reactions

Result = The pink material is obtained bakelite just.

Yield Procedure = Place 5 ml of glacial acetic acid and 2.5 ml of formaldehyde solution in a 500 ml beaker and add 2gm of Phenol. Then add conc. HCl dropwise into the mixture carefully. Keep the above mixture for 5 minutes. Firstly it will form white ppt later on pink colour hard mass is obtained. Wash the hard mass with distilled water to remove excess of impurities. Then remove the hard mass into the watch glass, dry and weight.

Result = Pink material is obtained


yield = 30 gm

Teacher's Signature _____

Object = To prepare urea formaldehyde (UF) resin

Reagent and apparatus = Urea formaldehyde, conc.
 H_2SO_4 , etc.

Reactions-



urea formaldehyde (UF) resin

Procedure = Place 20 ml of 40% formaldehyde soln in 500 ml Beaker. Add about 10gms of Urea while stirring until a saturated soln is obtained. Add few drops of conc. H_2SO_4 . Stim. continuously during the addition. All the

Teacher's Signature _____

voluminous white solid mass appears in the beaker. When rxn is complete, wash the residue with water and dry the product.

Result = white colour solid plastic.

yield = 28.18 gm

- Precautions = ① The glass apparatus should be cleaned and rinsed properly.
② H₂SO₄ should be mixed drop wise and stir continuously during the addition.
③ The reaction is sometimes vigorous and it better to be a few feet away from the beaker while adding the HCl and until the rxn is complete.

9 ~~10~~

Teacher's Signature _____

Experiment - 9

Object = To detect the elements and functional groups in given organic compound.

Detection of elements - Generally the organic compn have, C, H and oxygen. In addition to these three elements halogen, Nitrogen, and sulphur may also be present. Firstly these elements are converted or converted into ionisable in organic salts by fusion of the organic compound with metallic sodium. Detection of Nitrogen, sulphur, and halogens by this method is known as Lassaigne's test.


Apparatus

Chemicals Used = Test tube, ignition tube, tong and porcelain dish, beaker, funnel etc..


Chemicals used = Sodium metal, ferrous sulphate solⁿ, dil H₂SO₄, FeCl₃ solution, Sodium nitro nitrophthalide solⁿ, dilute acetic acid, lead acetate solⁿ, Cone, HNO₃, AgNO₃ solⁿ, dil HNO₃, chloroform, carbon tetrachloride, chlorine water, dil. HCl, ethyl ~~not~~ alcohol, acetone, cobalt nitrate solution.

Lassaigne's Test = Take a piece of dry sodium metal into a clean ignition tube and cover it with a piece of organic compound. Heat.

Teacher's Signature _____



($X = Cl, Br, F$)



If N and S are both present together




gently in the flame initially and then strongly until the bottom of the tube becomes red hot. Plunge the hot tube into a porcelain dish containing 15 ml of cold distilled water. Grind the contents with a pestle, boil the mixture for 5 mixture minutes and filter. The filtrate is known as sodium extract, which is then used for elemental detection.


Reactions Involved =

Reactions involved


Nitrogen:



ferrous hydroxide (dirty green)




Sodium ferricyanide



Ferricyanide (Prussian blue)

Both N and S:



ferric Sulphocyanide (blood red)





Sodium Nitroprusside

Sodium Eosulphenitroprusside (violet)



Teacher's Signature _____

Halogens:


(1) Chlorine -





(2) Bromine:



(3) Iodine:



(4) Iodine and Bromine =



Detection of functional group - In the identification of organic compound the detection of functional groups is an imp. step. The organic compounds

are divided into groups that exhibit similar chemical reactions and each group has a structural feature common to all organic molecules. This common structural feature is known as functional group.

Apparatus used = Beaker, test tube, boiling tube, funnel, burner, tripod stand, wire gauge, etc.

Chemicals used = Sodium bicarbonate solution, Natural ferric chloride solution, sodium nitrite, Sulphuric acid, Sodium hydroxide, phthalic anhydride, Sodium metal, Na_2SO_4 , Glacial acetic acid, ceric ammonium nitrate, ethanol, 2,4,-dinitrophenyl hydrazine, β -naphthal, CHCl_3 , alc. KOH , phenol, etc.

Table : 1

Sno	Experiment	Observation	Inference
	Carboxylic Groups ($-\text{COOH}$)		
i)	Litmus test = Place a crystal or solid sample or a drop of the given liquid sample on moist blue litmus paper.	Turns red	$-\text{COOH}$ group present.
ii)	Sodium bicarbonate test = take 2-3 ml of 10% solution of	Bubbles	$-\text{COOH}$ group present.

Teacher's Signature _____


sodium bicarbonate in a dry test tube add a pinch of given compound.

- (ii) Ester formation = Take small amount of sample compound in a dry test tube, add 2ml of ethyl alcohol and add 2-3 drops of conc. H_2SO_4 and warm on a water bath
- A fruity smell - $-COOH$ group present.


Reactions involved -:

(i) Litmus test: Blue \rightarrow Red

(ii) Sodium Bi-carbonate test:



(iii) Ester formation:



Phenolic Group -

- i) Neutral FeCl₃ test - Take 1 ml of neutral ferric chloride solution in cleaned and dry test tube and add 3-4 drops of a little bit of the organic compound. Shake well.
- Violet / red / blue / Green / Colour Phenolic Group present.


Teacher's Signature _____

ii)	litmus test = put a pinch of the organic compound on a moist blue litmus paper.	Turns red	Phenolic group present
iii)	Pthalic test = Take .5g of organic compound in a cleaned and dry test tube, add 1 gm of pthalic anhydride and 2 drops of concentrated H_2SO_4 . Heat gently for 1 minute. Cool and pour the mixture into a beaker containing 40ml of water and make it alkaline by adding NaOH.	Red / pink Resorcinol α -naphthol β -naphthol	phenol, o-cresol, catechol m-cresol Fluorescent Green, Green, Paint Green.
iv)	Liebermann's Nitroso test = Heat a small amount of the organic compound and few crystals of $NaNO_2$. Cool and add 1ml of concentrated H_2SO_4 and excess of water and excess of NaOH solution.	A deep green to blue soln is formed which turns red. Deep green to blue colour appears	

Reactions involved:

i) Litmus test: Blue \rightarrow Red

ii) Neutral FeCl₃ test \rightarrow






Teacher's Signature _____

Table 2

Table - 3

S.no	Experiment	Observation	Inference
	<u>Aldehyde group (-CHO)</u>		
i)	Schiff's Solution test = Take .5 gm of solid compound or 4 drops of liquid Compound in a test tube, add 1ml of Schiff's reagent, Shake vigorously.	Red or violet colour	-CHO Group present
ii)	Tollen's Reagent test = Take a little, Silver quantity of the organic comp ⁿ in mirror a test tube, add 2ml of Tollen's reagent, Shake warm and allow the ppt content to stand for 3-4 minutes	grey	-CHO Group present
iii)	Fehling's Solution test = Take 1ml of a mixture of equal amounts of Fehling soln A and B, add 4 drops or 0.2 gm of organic compound, boil for about 5 minutes	Red ppt	-CHO Group present

Reactions involved:

Teacher's Signature _____

Table - 5

S.no	Experiment	Observation	Inference
<u>Ketonic group -</u>			
i)	Schiff's Soln test = Take .5 g of solid comp ⁿ or 4 drops of liquid comp ⁿ in a test tube add 1ml of Schiff Reagent, shake vigorously.	No colour	Ketonic Group present.
ii)	2,4, Dinitrophenyl Hydrazine test = Take a little quantity of the organic comp ⁿ in a test tube add 2ml of 2,4-dinitrophenyl hydrazine reagent, shake, heat & cool.	A red, yellow or orange Coloured ppt	Ketonic Group present.

Table - 6

S.no	Experiment	Observation	Inference.
<u>Nitro Group (-NO₂)</u>			
i)	NaOH soln test - Take .5 gm of solid Comp ⁿ or 4 drops of liquid comp ⁿ in a test tube add 1ml of ag. soln of NaOH and warm.	Yellow Colour	Nitro Group present.
ii)	Mulliken's test - Take a little quantity of org. Comp ⁿ in the dry test tube add 1ml of ethanol, 1ml of CaCl ₂ or NH ₄ Cl soln & a pinch of zinc dust, boil over content for 5 min. Cool & filter into a test tube contain 2ml of toluidine reagent.	Grey or Black ppt	Nitro Group present.

Teacher's Signature _____

Table - 7

Sno	Experiment	Observation	Inference
①	Primary amines (-NH ₂)		
i)	Litmus test - Put a pinch of the organic comp ⁿ on a moist red litmus paper.	Turns blue	Amines Group present.
ii)	Nitrous acid test = Dissolve a small Brisk quantity of organic comp ⁿ in effervescence dil. HCl and cool. Add 10% ag. NaNO ₂ Solution.	Aliphatic amines primary amines Group present	
iii)	Diazotisation test - Dissolve about 2 gm of the organic comp ⁿ in dil HCl and cool. Add 10% ag. NaNO ₂ soln. Pour the reaction content in a Beaker containing alkaline β -naphthol solution	Red or orange dye	Aromatic primary amines Groups present

Teacher's Signature _____

Experiment - 10

Object Result = we have detect the following functional groups in the given organic comp'

- ① Carboxylic group (-COOH)
- ② Alcoholic group (-OH)
- ③ Phenolic group
- ④ Aldehydic Group (-CHO)
- ⑤ ketonic Group present
- ⑥ Nitro group (-NO₂)
- ⑦ Primary amines (-NH₂)

Precautions =

- ① Acid burns = Wash with plenty of water then with very dilute Sodium bisulphite solution
- ② Bases = Wash with plenty of water then with very dilute acetic acid soln.
When you are heating a test tube, never point its mouth towards yourself or anyone else.
- ③ Protect your hands if you have to pick up hot objects.
- ④ Wash test tube before adding.

⑨ Praj

Teacher's Signature _____

Experiment - 10


Object = To determine the strength of given hydrochloric acid solution by titrating it pH metrically against Sodium Hydroxide.

Chemicals / Apparatus Required = pH meter, glass electrode, reference electrode, burette, stirrer, NaOH, HCl.

Principle = When an alkali is added to an acid solution the pH of the solution increases slowly. But about the equivalence point, the rate of change of pH of the solution is very rapid. From the sharp break in the curve, we can find the equivalence point, from which the strength can be calculated by Normality equation.

Procedure = First standardize the PH meter against a buffer soln of known pH. Now first wash the glass electrode and reference electrode with distilled water and then rinse with the acid solution. Take 50ml of HCl solution in a 500 ml of beaker. Add sufficient water so that the glass electrode as well as reference electrode is completely dipped. Note the pH of pure acid soln. Now add 1ml of 1N NaOH from the burette in the beaker. Stir the

Teacher's Signature _____



the contents well and note the pH of the solution. Continue adding NaOH soln from the Burette and note down the pH after each addition. Near the equivalence point the change in pH occurs rapidly then add alkali in fraction.

Observations = Volume of acid taken = 50 ml (say)

Volume of NaOH added (ml)	0	1	2	3	4	5	6	7	8	9	10	11	12	13
pH	0.60	0.62	0.40	0.68	0.72	0.78	0.86	0.98	1.14	1.60	9.14	9.80	9.92	

Calculations = A curve is plotted with pH values as ordinate and volume of NaOH added ml as abscissa. The sharp break in the curve corresponds to the equivalence points. Suppose Vol. of 0.1 NaOH used at end point = V₂ ml
(Acid) (alkali)

$$N_1 V_1 = N_2 V_2$$

$$N_1 \times 50 = \frac{1}{10} \times V_2$$

$$N_1 = \frac{\frac{1}{10} \times V_2}{50} = \frac{6.5}{500} = 0.013$$

$$\text{Strength} = N \times \text{Eq. wt.}$$

$$= 0.013 \times 36.5 \text{ gm/l} = 0.4795 \text{ gm/lit.}$$

Teacher's Signature _____

Result = The strength of given HCl soln is
• 4-45 gml lit

Precautions = The temperature of central knob of the pH meter should be adjusted to the room temperature.

- After the addition of alkali, the soln should be thoroughly stirred
- The pH meter should be first standardized by taking a buffer of known pH.

Viva Voce

Q1- Define pH?

Ans- pH is the scale to specify how acidic or basic a water-soln is.

Q2- Do you use any indicator in pH titration - No.

Q3- What is the desirable pH range for drinking water?

Ans- 6-7.5 desirable range.

Q4- What is effect of temperature on pH?

Ans- pH decreases with increase in temperature

Q5- Who suggested the pH value

Ans- Danish chemist Soren Pedar Lauritz Sorenson

(B) 

Teacher's Signature _____